



# HeartSpan<sup>®</sup>

Fixed Curve Braided Transseptal Sheath



# HeartSpan®

Fixed Curve Braided Transseptal Sheath

## Product Features

Braided steel core imparts excellent torqueability

Lubricious coating inside and out

Dual distal holes

Tapered dilator and tip geometry enhance tissue crossing

Snap-fit dilator locks in place

Radiopaque sheath and dilator enable in situ visualization

Soft, atraumatic tip

Embedded platinum marker

Works in unison with the HeartSpan Transseptal Needle

HeartSpan Fixed Curve Kit



- (1) Braided Transseptal
- (1) Dilator
- (1) J-tipped Guide wire

### HEARTSPAN Fixed Curve Transseptal Sheaths

Part Number	French Size	Length	Curve
FCL-160-00	8.5F	60 cm	15° Curve
FCL-160-01	8.5F	60 cm	30° Curve
FCL-160-02	8.5F	60 cm	55° Curve
FCL-160-03	8.5F	60 cm	120° Curve (short)
FCL-160-04	8.5F	60 cm	120° Curve (long)
FCL-160-05	8.5F	60 cm	150° Curve
FCL-160-06	8.5F	79.4 cm	15° Curve
FCL-160-07	8.5F	79.4 cm	55° Curve
FCL-160-08	8.5F	79.4 cm	90° Curve
FCL-160-09	8.5F	79.4 cm	120° Curve (long)
FCL-160-10	8.5F	101.5 cm	55° Curve

