

Worley™ VEIN SELECTOR



AERTMEDICAL™

Advancing Vessel Selection

Product Features

Unique accordion hub design to minimize kinking and enhance steerability

Proprietary braided shaft design to maximize flexibility

Soft tapered tip to reduce vessel trauma


Highly visible under fluoroscopy

Improve LV lead location.¹

¹ Jackson KP, Hegland DD, Frazier-Mills C, Piccini JP, Koontz JI, Atwater BD, Daubert JP, Worley SJ. Impact of using a telescoping-support catheter system for left ventricular lead placement on implant success and procedure time of cardiac resynchronization therapy. Pacing Clin Electrophysiol. 2013 May; 36(5):553-8.

WORLEY VEIN SELECTOR

Part Number	French Size	I.D.	Length	Curve
57538CS-WOR	5F	.046"/ 1.17 mm	75 cm	Standard
57538CSV-WOR	5F	.046"/ 1.17 mm	75 cm	Vertebral
57538CSHK-WOR	5F	.046"/ 1.17 mm	75 cm	Hook



Greater visibility for reliable access to the target vein



Ergonomic wing hub enhances handling and control

